



# Microbiome pipelines : 16S sequencing pipeline


[jianhong/16s\\_pipeline](#)




# stats




# Metagenomics




# Microbiome profiling



# Microbiome








- Template relatively short
  - 1.5 kb
- Highly conserved and well studied
  - Amplified by same primers
- Genus-level resolution
  - Closely-related species have a high sequence similarity across the 16S gene

# 16S rRNA-seq

515F primer



806R primer



primers

4 / 27 / 22

Primer name	Primer sequence	Reference
515f Original	GTGCCAGCMGCCGCGGTAA	Caporaso et al.
806r Original	GGACTACHVGGGTWTCTAAT	Caporaso et al.
515f Modified	GTGYCAGCMGCCGCGGTAA	Parada et al.
806r Modified	GGACTACNVGGGTWTCTAAT	Apprill et al.
926r	CCGYCAATTYMTTTRAGTTT	Parada et al.
ITS1f	CTTGGTCATTTAGAGGAAGTAA	Gardes and Bruns
ITS2	GCTGCGTTCTTCATCGATGC	White et al.

HTTPS://WWW.IDT.DNA.COM/PAGES/EDUCATION/DECODED/ARTICLE/16S-RNA-INDEXED-PRIMERS-AMPLIFY-PHYLOGENIC-MARKERS-FOR-MICROBIOME-SEQUENCING-ANALYSIS

# Note

- There are multiple different pipelines, such as mothur, QIIME, DADA2, and etc.
- Basic steps: sequence trimming -> demultiplexing -> Chimera filtering -> Sequence classification
- There are multiple available 16S databases, and the used database will affect the taxonomic classification.



# Current pipeline

BCL2fastq

- bcl2fastq

Trimming

- Trimmomatic

demultiplex


- fastq\_pair\_filter.py
- QIIME2::demux

Filter reads

- DADA2

classification

- DADA2



# nextflow

- Track progress tracking
  - Troubleshooting
  - Version control
- Module based design
- Portable and scalable
- Highly reproducible
- Compatible with cluster/cloud computation

# Why CONDA®

- Conda is an open-source package management system and environment management system
- It does not require root privilege
- Cons:
  - Not stable as Docker/Singularity
  - Version conflicts
  - Storage space

```
(base) ouj@Jianhongs-MacBook-Pro-2 tmp4genomictools % which ssh
/usr/bin/ssh
(base) ouj@Jianhongs-MacBook-Pro-2 tmp4genomictools %
```

## 3 Steps to Run jianhong/16S\_pipeline

- Install `conda`
- Install `nextflow`
- Run pipeline



<https://youtu.be/XiVnM5iptUI>

# Install `MiniConda`

- `wget -O miniconda.sh "https://repo.anaconda.com/miniconda/Miniconda3-latest-Linux-x86_64.sh"`
- `bash miniconda.sh`
- `vim ~/.condarc` (The order is important)
  - channels:
    - conda-forge
    - bioconda
    - defaults

# Create `nextflow` environment

- `conda create -y --name nextflow bioconda::nextflow=21.10.6`

# Test run

- Create an interactive job in DCC
  - `srun --mem 10G -c 2 --pty bash -i`
  - `conda activate nextflow`
  - `nextflow run jianhong/16S_pipeline -r main -profile conda,test`
  - To save time for testing:
    - `--silva_nr99 '/work/jo117/16S_pipeline/silva_nr99_v138.1_train_set.fa.gz'`
    - `--silva_tax '/work/jo117/16S_pipeline/silva_species_assignment_v138.1.fa.gz'`

# Test run by submit it as a job

- Create profile config file named as profile.config

```
// submit by slurm
process.executor = "slurm"
process.clusterOptions = "-J ProjectName"
params {
  max_cpus = 2
  max_memory = '6.GB'

  // Input data
  input = "${projectDir}/assets/test_data"
  skip_bcl2fastq = true
  barcodes = "${projectDir}/assets/barcodes.tsv"
  metadata = "${projectDir}/assets/metadata.csv"

  // report email
  email = 'your@email.addr'
}
```

- Create a slurm script file named as microbiome.sh

```
#!/bin/bash
#SBATCH -J 16S_submitter #jobname
#SBATCH -o microbiome.out.%A_%a.txt
#SBATCH -e microbiome.err.%A_%a.txt
#SBATCH --mem-per-cpu=10G #memory for the job submission
node
#SBATCH -c 1 # 1 CPU is good enough

source ${HOME}/.bashrc
conda activate nextflow

nextflow run jianhong/16S_pipeline -r main -profile conda -c
profile.config -resume
```

**sbatch microbiome.sh**



# Run pipeline for your own data

- Create profile config file named as profile.config

```
// submit by slurm
process.executor = "slurm"
process.clusterOptions = "-J ProjectName"
params {
// Input data
input = 'path/to/your/initialFiles' // replace it by your own folder contain
Intensities folder.
barcodes = 'path/to/your/barcodes.tsv'
metadata = 'path/to/your/metadata.csv'
// report email
email = 'your@email.addr'
}
```

- Create a slurm script file named as microbiome.sh

```
#!/bin/bash
#SBATCH -J 16S_submitter #jobname
#SBATCH -o microbiome.out.%A_%a.txt
#SBATCH -e microbiome.err.%A_%a.txt
#SBATCH --mem-per-cpu=20G #memory for the job submission node
#SBATCH -c 1 # 1 CPU is good enough

mkdir -p tmp
export TMPDIR=${PWD}/tmp
export TMP=${PWD}/tmp
export TEMP=${PWD}/tmp
source ${HOME}/.bashrc
conda activate nextflow
module load bcl2fastq/2.20

nextflow run jianhong/16S_pipeline -r main -profile conda -c profile.config
-resume
```

**sbatch microbiome.sh**

# Parameters: input/output


- `--input '[path to raw reads files]'`
  - If fastq files are supplied, the names must follow the pattern: `_[RI][12]_[0-9]+`, eg: `sample1_R1_001.fastq.gz`
- `--barcodes '[path to barcodes tsv file]'`
  - Barcodes will be used by QIIME2::demux
  - Sample file can be found at [https://github.com/jianhong/16S\\_pipeline/blob/main/assets/barcodes.tsv](https://github.com/jianhong/16S_pipeline/blob/main/assets/barcodes.tsv)
- `--metadata '[path to metadata csv file]'`
  - Metadata will be used by Bioconductor::phyloseq package, see help <https://bioconductor.org/packages/phyloseq/>
  - Sample file can be found at [https://github.com/jianhong/16S\\_pipeline/blob/main/assets/metadata.csv](https://github.com/jianhong/16S_pipeline/blob/main/assets/metadata.csv)

# Parameters: 16S reference

- ❑ --silva\_nr99: used by assignTaxonomy
- ❑ --silva\_tax: used by addSpecies
- The latest release can be found at: <https://www.arb-silva.de/download/arb-files/>
- The defaults are set as files in <https://zenodo.org/record/4587955#.YiZhRpNudJU>
- The references will be used by DADA2::assignTaxonomy and addSpecies, see help at <https://benjjneb.github.io/dada2/training.html>
- They can be replaced by other resources such as RDP and UNITE.

# Parameters: pipeline control

- `--bcl2fastq`: use DCC module `bcl2fastq/2.20`
- `--skip_bcl2fastq`: skip `bcl2fastq` or not
- `--skip_fastqc`: skip `fastqc` or not
- `--skip_demultiplex`: skip `demultiplex` or not



# Errors

- Exit code: 127: command not found
  - Check the conda environment installation
- Exit code: 137: out of memory
  - Increase the resource requirement
- Exit code: 139: fault installation
  - Check the conda environment installation

# Increase the resource requirement

- [https://jianhong.github.io/16S\\_pipeline/usage.html#resource-requests](https://jianhong.github.io/16S_pipeline/usage.html#resource-requests)
- The maximal limitation
  - max\_cpus, max\_memory, max\_time in profile.config file, eg:  
[https://github.com/jianhong/16S\\_pipeline/blob/main/conf/test.config](https://github.com/jianhong/16S_pipeline/blob/main/conf/test.config)
- Change the requirement by process name in profile.config file.
  - `process { withName: BCL2FASTQ { memory = 100.GB } }`

# Module specific parameters

- All default settings can be found at [https://github.com/jianhong/16S\\_pipeline/blob/main/conf/modules.config](https://github.com/jianhong/16S_pipeline/blob/main/conf/modules.config)
- Eg: for trimmomatic, see help at [http://www.usadellab.org/cms/uploads/supplementary/Trimmomatic/TrimmomaticManual\\_V0.32.pdf](http://www.usadellab.org/cms/uploads/supplementary/Trimmomatic/TrimmomaticManual_V0.32.pdf)

```
withName: REMOVE_PRIMERS {  
    ext.args = 'ILLUMINACLIP:techseqs.fa:2:30:10'  
    publishDir = [  
        path: {"${params.outdir}/1_remove_primers"},  
        mode: 'copy',  
        saveAs: { filename -> filename.equals('versions.yml') ? null : filename }  
    ]  
}
```

# Parameters with demuxing

```
withName: QIIME_DEMUX {
```

```
  ext.args = '--m-barcodes-column barcode-sequence --p-rev-comp-mapping-barcodes '
```

```
  publishDir = [
```

```
    path: { "${params.outdir}/3_demultiplex" },
```

```
    mode: 'copy',
```

```
    saveAs: { filename -> filename.equals('versions.yml') ? null : filename }
```

```
  ]
```

```
}
```

sample-id	barcode-sequence
#q2:types	categorical
Spinach1	TGTGCGATAACA
Spinach2	GATTATCGACGA
Spinach3	GCCTAGCCCAAT

<https://docs.qiime2.org/2022.2/plugins/available/demux/emp-paired/>



# Trim by dada2::filterAndTrim

For `FILTERING`, the options are

<code>--trimming_reads, -t</code>	<code>"logical",</code>	Trim reads or not.
<code>--trim_left, -a,</code>	<code>"integer",</code>	Default 0. The number of nucleotides to remove from the start of the R1 read. If both <code>trunc_length_left</code> and <code>trim_left</code> are provided, filtered reads will have length <code>trunc_length_left-trim_left</code> .
<code>--trim_right, -b,</code>	<code>"integer",</code>	Default 0. The number of nucleotides to remove from the start of the R2 reads. If both <code>trunc_length_right</code> and <code>trim_right</code> are provided, filtered reads will have length <code>trunc_length_right-trim_right</code> .
<code>--trunc_length_left, -m,</code>	<code>"integer",</code>	Default 0 (no truncation). Truncate R1 reads after <code>trunc_length_left</code> bases. Reads shorter than this are discarded.
<code>--trunc_length_right, -n,</code>	<code>"integer",</code>	Default 0 (no truncation). Truncate R2 reads after <code>trunc_length_right</code> bases. Reads shorter than this are discarded.

# dada2 workflow

SOURCE:  
HTTP://BENJNEB.GITHUB.IO/DADAD2\_PIPELI  
NE\_MV/BIGDATA.HTML

1. check PAIRED END by `length(filtered forward) == length(filtered reverse)`
2. Learn Error Rates by a subset data (36 samples) (Another way `learnErrors`)
  1. Dereplication by `dada2::derepFastq`
  2. Sample Inference by `dada2::dada(dereps, err=NULL, selfConsist=TRUE, multithread=NCORE)`
3. Run `dada` by the error model for all samples
  1. `derepFastq` → `dada(sample, err=errModel, multithread=NCORE)`
4. Merge Paired Reads by `mergePairs`
5. Construct Sequence Table: `makeSequenceTable` and Trim sequences by sequence lengths.
6. Remove Chimeras: `removeBimeraDenovo(seqtab, method='consensus', multithread=NCORE)`
7. Assign Taxonomy:
  1. `TaxTable <- assignTaxonomy( sequences, TRAIN_SET, tryRC=TRYRC, multithread=NCORE)`
  2. `addSpecies(TaxTable, SPECIES_ASSIGNMENT, tryRC=TRYRC, verbose=TRUE)`

# Parameter - - tryRC


- If you have unmerged reverse-complement reads issue, you may want to try the parameter --tryRC
- It will try to merge the reverse-complement reads counts.
- Please note that partial of this operation is not in DADA2 package. It is used for the reads that not proper demultiplexed.

# Modula names and outputs

Modules	program	outputs
BCL2FASTQ	bcl2fastq	0_data_raw/*.fastq.gz
FASTQC	fastqc	fastqc/*_fastqc.{html,zip}
REMOVE_PRIMERS	trimmomatic	1_remove_primers/*.fastq.gz
SYNC_BARCODES	sync_paired_end_reads.py	2_sync_barcodes/*.fastq.gz
QIIME_DEMUX	qiime demux	3_demultiplex/demuxd_reads/*.fastq.gz
FILTERING	dada2::filterAndTrim	4_filter/<sampleid>/*.fastq.gz
DADA2	dada2	5_data2/*
FHYLOSEQ	phyloseq	6_phyloseq/*
KRONA	Krona	7_Krona/*
MULTIQC	multiqc	multiqc/*


# Krona visualization

- [https://jianhong.github.io/16S\\_pipeline/krona.html](https://jianhong.github.io/16S_pipeline/krona.html)
- <https://github.com/marbl/Krona/wiki>





# QC reports

- [https://jianhong.github.io/16S\\_pipeline/multiqc.html](https://jianhong.github.io/16S_pipeline/multiqc.html)



# Get help

- Online Doc: [https://jianhong.github.io/16S\\_pipeline/](https://jianhong.github.io/16S_pipeline/)
- Email: [jianhong.ou@duke.edu](mailto:jianhong.ou@duke.edu)
- Report an issue: [https://github.com/jianhong/16S\\_pipeline/issues](https://github.com/jianhong/16S_pipeline/issues)
- Dcc usage: <https://oit-rc.pages.oit.duke.edu/rcsupportdocs/dcc/#cluster-shared-storage-resources-work-and-scratch>



## Acknowledgements

